
truight

Trylight for your biosensing!





Andra da in

\Diamond

MAIN BIOSENSING APPLICATION

- Magnetoelastic sensor
- for reproducibility
- Quartz Crystal Microbalance (QCM) *for high specificity*
- Surface Plasmon
 Resonance (SPR)
 for kinetics studies
- Plasmonic Enhanced Fluorescence (PEF) for high sensitivity

The Photochemical Immobilization Technique (PIT) is a functionalization method able to tether antibodies (Abs) on metal surfaces (such as gold or other bio-compatible metals) in a proper orientation that is with their binding sites well exposed to the environment.

This strategy is based on the selective photochemical reduction of disulfide bridges in immunoglobulins (Ig) produced by UV activation of near aromatic amino acid. Briefly, PIT is based on the selective photoreduction of disulphide bridges produced by the UV activation of the trp/cys-cys by means of a Trylight lamp.

The UV photon energy is adsorbed by tryptophan and transferred to near electrophilic species like the close cys-cys resulting in the cleavage of the disulphide bridge and in the formation of new reduced thiol (SH) groups able to bind. Considering that the triad of residues trp/cys-cys can be found in every Ig, the PIT is applicable in a wide range of fields.

How Trylight works


1 A solution of IgG

A solution of IgG antibodies is pipetted in a quartz cuvette

•••••••


2

The quartz cuvette is placed inside the device



The UV irradiation of 30 s is started

2



lgG

The irradiated solution can be convoyed onto any gold transducer by fluidics and/or dropping

lgG

4

$\sqrt{}$

Technical sheet

LAMP TYPE PR7GL-5 253.7 nm

LAMP BODY Pure quartz of European origin

PIPE DIAMETER 10 mm

FORM U

LENGTH total pin body excluding 45 mm

WIDTH 40 mm

SUPPLY VOLTAGE 12 Vdc

LAMP CURRENT 170 mA

POWER OF SINGLE LAMP 6 W

WORKING FREQUENCY 35 KHz

UV EMISSION WAVELENGTH TYPE PR7GL-5 254 nm

TOTAL UV INTENSITY ON THE CUVETTE 0.3 W / cm²

OPERATING TEMPERATURE -10 to 70 °C

HUMIDITY max 95%

CONNECTION FEET tungsten diam. 0.5 m

trylight

How to order You can order your Trylight directly online at <u>www.trylight.eu</u> or by e-mail at <u>info@trylight.eu</u>

Cost The cost of Trylight is 1400 Euro, vat excluded

Payment

You can pay online via paypal with most of credit cards or you can pay by bank transfer. Ask further details at info@trylight.eu

Shipping

We can ship Trylight worldwide. Ask for a specific quotation at info@trylight.eu

www.trylight.eu info@trylight.eu

Trylight is patented. Utility Model Patent nr. 202018000003368 Trylight is a product of



Promete srl Piazzale Tecchio 45 80125 Napoli tel +39 081 056850/51 promete@promete.it

Realized under the project NEST that has received funding from the Regione Campania POR CAMPANIA FESR 2014-2020 ASSE III, OBIETTIVO SPECIFICO 3.4. AZIONE 3.4.2. Avviso pubblico per la concessione di contributi finalizzati al finanziamento di programmi di internazionalizzazione delle Micro e PMI campane.







